DEVELOPMENT OF AN ELISA FOR THE DIRECT MEASUREMENT OF FREE 250H VITAMIN D

N. Heureux^{*}, M. Anciaux, M. Poncelet, F. Mathieu DIAsource Immunoassays, Louvain-la-Neuve, Belgium L.M. Swinkels, T. Huijs, S. Boerrigter, E. Lindhout, G. Mayer, M. Martens Future Diagnostics B.V., Wijchen, The Netherlands

1. Introduction

Almost all circulating 250H Vitamin D in serum is bound to Vitamin D Binding Protein DBP (88%) and Albumin (12%). A very small fraction, approximately 0.04% of the 250H Vitamin D, circulates in the free, non-protein bound form.¹

The concentration of DBP is not constant and can be influenced by a number of factors including obesity,² pregnancy,³ the use of oral contraceptives,⁴ hormone replacement therapy,⁵ liver disease,⁶ renal disease,⁷ proteinuria⁸ and intensive care.⁹

Furthermore, common genetic polymorphisms in the DBP gene produce variant proteins that differ

4. Assay Performance

Range	0-40 pg/mL
Precision	<10%
LoB	0.7 pg/mL
LoD	1.9 pg/mL

Additional experiments have shown that the addition of albumin or Vitamin D Binding Protein to

in their affinity for 250H Vitamin D.^{10,11}


For thyroid and steroid hormones it has been shown that their *in vivo* physiological activity correlates better with the free, non-protein bound fraction, than with the total concentration of the hormone in plasma. In such conditions the measurement of the concentration of the free circulating hormone provides more physiologically relevant information. This notion is known as the "free hormone hypothesis".^{12,13} Particularly in situations in which the level of binding proteins is elevated or decreased, the measurement of total circulating hormone may lead to a wrong diagnosis.

The bioavailability of 250H Vitamin D may be more complex, since it is a precursor that is metabolized in the kidney to the active 1,25(OH), Vitamin D. However, literature data suggest a role for the measurement of free 25 OH Vitamin D.¹⁴⁻¹⁷

Here we describe an ELISA that offers the possibility to measure the concentration of free 250H Vitamin D in serum.

2. Principle of the Assay


The Free 250H Vitamin D ELISA is based on a two-step immunoassay procedure. During the first incubation step free 250H Vitamin D binds to the monoclonal anti-Vitamin D in the microtiter plate. The in *vivo* equilibrium between free and bound 250H Vitamin D is minimally disturbed by the use of a specific displacement reagent. After washing, a fixed amount of biotinylated 250H Vitamin D is added to each well. The non-bound biotinylated 250H Vitamin D is removed by washing and a streptavidin peroxidase conjugate



serum leads to a decrease in the observed level of free 250H Vitamin D.


Sample	Free 250H Vitamin D (pg/mL)
Native sample 1	3.25
Native sample 1 + 60 g/L Albumin	1.95
Native sample 2	13.0 pg/mL
Native sample 2 + 1 g/L DBP	5.8 pg/mL

5. Pregnancy and Renal Failure



Box and Whisker Plot of Total 25OH Vitamin D concentration by LC-MS

-----95% CI Notched Outlier




Box and Whisker Plot of DBP concentration by ELISA


Total 250H Vitamin D, DBP and free 250H Vitamin

is added. In a next step TMB chromogenic substrate is added. The reaction is then stopped by addition of Stop solution and the absorbance is measured using a plate reader.

The concentration of free 250H Vitamin D in the calibrators was determined using a rate dialysis technique.

3. Assay Protocol





Box and Whisker Plot of free 250H Vitamin D concentration by **ELISA**

6. Discussion


concentrations were measured in normal controls "n" (n=45), Pregnant "p" (n=31) and Renal Failure "r" patients (n=9).

250H Vitamin D levels are higher in Total pregnancy than in controls. However, levels of free 250H Vitamin D suggest that the Vitamin D status in pregnancy is lower than in controls. The higher concentration of DBP in pregnancy samples corroborates these results.



Join the Scientific Community who talks about free Vitamin D: LinkedIn group "Free [™] Vitamin D″









DOI: 10.3252/pso.eu.17ece.2015